

SAPIENS-1MS: Ultra Inert Test (TK-SAP1MS) Performance against major ultra inert column manufacturers

Column: SAPIENS-1MS P/N TR-430232
Dimensions: 30m x 0.25mm x 0.25µm
Injection: 0.1µL TEST ULTRAINERT, split 500:1
Injector temperature: 275°C
Detector temperature: 340°C
Carrier gas: H₂ 12psi (82.8kPa)
Oven temperature: 65°C (Isothermal)
Detector: FID, 340°C

