

**SAPIENS-5MS: Ultra Inert Test (TK-SAP5MS)**
**Performance against major ultra inert column manufacturers**

Column: SAPIENS-5MS P/N TR-420232


Dimensions: 30 m x 0,25 mm x 0,25µm

Injection: 0,1µL TK-SAPX5MS, Split 500:1, 275°C

 Carrier gas: H<sub>2</sub>, 12 psi (82,8 kPa)

Oven temperature: 65°C (isothermal)

Detector: FID, 340°C


**SAPIENS-5MS: Ultra Inert Test (TK-SAP5MS)**
[www.teknokroma.com](http://www.teknokroma.com)