

SAPIENS-5MS: Acidity Test - Perfect peak shapes

Column: SAPIENS-5MS P/N TR-420232


Dimensions: 30 m x 0,25 mm x 0,25µm

Injection: 0,3µL Phenols EPA 604 (2 mg/ml MeOH), Split 100:1, 275°C

Carrier gas: H₂ Ct flow 1,4 ml/min

Oven temperature: 80°C (4 min) to 200 °C (2 min) @ 8°C/min

Detector: FID, 300°C



SAPIENS-5MS: Acidity Test

www.teknokroma.com