

## SAPIENS-WAX.MS: Ylang-Ylang Oil by GC-FID

Column: SAPIENS-WAX.MS P/N TR-480262


Dimensions: 60 m x 0,25 mm x 0,25 µm

Injection: Ylang-Ylang oil (neat wet needle), 275°C

Carrier gas: He, 25 psi (172.2 kPa)

Oven temperature: 40°C (hold 3 min) to 260°C (hold 30 min) @ 3°C/min

Detector: FID @ 295°C



SAPIENS-WAX.MS

[www.teknokroma.es](http://www.teknokroma.es)

Teknokroma Analitica S.A.- Camí de Can Calders, 14 - 08173 Sant Cugat del Vallès (Barcelona) Spain

Tel: +34) 936 748 800 Fax: +34) 936 752 405 e-mail: comercial@teknokroma.es